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Miplate provides optimized anatomical screw position, soft tissue problems reduction and there is no need for plate contouring.

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Angular stability

Miplate Prevents screw loosening as well as primary and secondary loss of reduction, preserves bone vascularization, allows early functional mobilization and shows a great improvement for osteoporotic bones.

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Less invasive procedure

A radiolucent handle ref.7027 facilitates the insertion of the plate as well as accurate and hassle-free percutaneous placement of the screws. Additional instrumentation facilitate indirect reduction.

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Anatomic reduction

Fixation of extra- and intra-articular distal femur fractures with the precontoured LCP DF plates allows for anatomic reduction.

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Stable fixation

Locking holes allow fixation with locking screws for angular stability. A fixed-angle construct is advantageous in osteoporotic bone and multifragment fractures where traditional screw purchase is compromised.

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Preservation of blood supply

The Miplate technique preserves the blood supply through a minimally invasive surgical technique and minimal bone-to-plate contact.

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Early, active mobilization

Miplate provides stable fracture fixation with minimal trauma to the vascular supply. This helps improve the environment for bone healing, accelerating the patient's return to previous mobility and function.

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


Use the pre operative template to determine the length of the plate and the screws positioning. Generally, the length of the plate and the screw positioning are chosen in the same way as an external fixator is. The template is also used to determine in which hole must the screw be inserted in the fracture region.

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
Preoperative screw-length selection using an AP radiograph

To select the proper screw length for the condyle, it is possible to perform a preoperative x-ray with the 50 mm wide calibrator and to use the table below.

1. Place the x-ray calibrator medially or laterally at the height of the condyle.
2. Take an AP radiograph of the distal femur.
3. Measure the width of the x-ray calibrator (XRC) on the radiograph.
4. Measure the maximum condyle width (MCW) on the radiograph.
5. Determine the real condyle width (RCW).
6. Determine the screw lengths for the screw holes A to G using the table below. The positions A to G are indicated on the x-ray template and on the Miplate insertion guide.

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Patient's Position

Position the patient supine on a radiolucent table. The leg should be freely movable. The contralateral leg can be placed in an obstetric leg holder. Place the knee joint line slightly distal to the hinged part of the table to allow flexion of the knee during surgery.

Avoid too strong a traction and a fully extended knee, as the forces of the gastrocnemius muscle would draw the distal fragment into recurvatum. This does not only make the reduction of the fracture difficult, but also endangers the popliteal artery and vein.

In very short distal fragments, it is recommended to flex the lower leg to approximately 60°. This also reduces the traction force of the gastrocnemius muscle.

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
Instruments Assembly

Join the two parts of the insertion guide (Pic.01 and Pic. 02). Insert the Cannulated guide Ø 7,5mm (short) Ref.7013 (pic.03) into the hole A in the insertion guide. Insert the guide with thread Ø 2,5 mm Ref.7008 into the guide Ref. 7013 (Pic.04).

Put the insertion guide in the three parts of the locking plate mechanism (pic.05). Thread the guide Ø 2.5mm Ref. 7008 into the plate tightening it in a smooth way with a curved wrench Ref.7009 (Pic.06) and traction the plate with the wrench Ref.7048 (Pic.07). For a better plate fixation in the insertion guide, introduce the guide with thread Ref.7007 or Ref.7008 and the Cannulated guide Ø 7,5 long Ref.7018 in the hole B and thread the plate.

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Note: To facilitate implant removal, close the unoccupied screw holes with screw hole inserts prior to inserting the plate. Use the torque-limiting screwdriver. The optimum torque is reached after one click.

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Reduce the fracture

In an intra-articular fracture, first reconstruct and stabilize the entire joint. The figure shows the possible positioning sites for lag screws in the condyles. Take care to ensure that these lag screws will not collide with the screws inserted through the insertion guide.

The fracture can be aligned manually by traction using a temporary knee-bridging external fixator or a distractor. Intraoperative x-ray or image-intensifier control is recommended to check the reduction. The anteromedial insertion of a Schanz screw can be beneficial in distal fragment manipulation.

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Extra-articular Fractures

Perform a skin incision of 8cm in the Gerdy's tubercle. Perform a skin incision from Gerdy's tubercle about 8 cm in a proximal direction. Open the space between the lateral vastus and the periosteum. The plate can be inserted into the space between the periosteum and the muscle.

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Intra-articular fractures

In intra-articular fractures, an anterolateral arthrotomy providing good control of the reduction is recommended. This arthrotomy also allows a subsequent insertion of the plate and can be used to insert lag screws from medially.

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Use the assembled insertion guide to insert the plate between the lateral vastus muscle and the periosteum. Slide the plate proximally and ensure that its proximal end remains in constant contact with the bone. Position the distal end of the plate against the lateral condyle. To find the correct position, move the plate proximally and then back distally until the plate fits the condyle (Pic.09).

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Due to its weight, the insertion guide tends to tilt dorsally. If the insertion guide points parallel to the floor with the patient in a supine position, the plate is externally rotated and no longer lies flat up against the lateral condyle. The fixation bolt must be oriented parallel to the patello-femoral joint. Consequently, the insertion guide shows an internal rotation of about 10°. This occurrence is also visible in the AP view of an image intensifier. The plate must lie flat up against the condyle to ensure an optimal fit on the bone.


Fig.09

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Once the plate is properly aligned with the bone, remove the drill sleeve and stabilization bolt from hole B -Ref.7007 or Ref.7008. Insert the trocar Ref.7012 through the long drill sleeve Ø 7,5 Ref.7018 in the most proximal hole of the plate (Pic.10). Perform a stab incision and push the drill sleeve and the trocar down to the plate (Pic.11). Check the correct position of the proximal part of the plate, either with the image intensifier or by direct palpation.


Fig.10


Fig.11

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Note: Due to soft tissues around the stabilization bolt Ref.7016, it will be difficult to change the position of the plate/handle assembly once the bolt has been inserted.

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Optional: Use a Kirschner wire to check the correct position of the proximal part of the plate on the bone.

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For preliminary fixation of the plate, use 2.0 mm Kirschner wires through the fixation and stabilization bolts. Carefully check the position of the plate and the length of the reduced injured limb. Once the reduction has been successfully completed and the plate has been positioned correctly, the locking screws can be inserted (Pic. 12).


Fig.12


Fig.13

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If necessary, it is possible to use 2.0 mm Kirschner wires for the preliminary fixation along the full length of the plate. Use the aiming device for Kirschner wires to insert the wires on the ventral and dorsal side of the plate. Note that the distance between bone and plate should be kept as short as possible when inserting the wires, as they are arranged in a convergent way. After the insertion of the Kirschner wires, the distance between plate and bone can no longer be reduced. After removing the Kirschner wire sleeves and the aiming device, proximal/distal displacement and adjustment of the position of the plate can be carried out. At the same time, the lateral Kirschner wires prevent the plate from migrating into the sagittal plane. Once the correct position is determined, the plate can be locked temporarily with a Kirschner wire through the fixation bolt (Pic. 13)

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Screw placement depends on the type of fracture. The position of the screws should be chosen in accordance with established biomechanical principles for internal fixation. The screws should be inserted close to and remote from the fracture gap in the main fragments. Use at least four screws per fracture side.

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Once the initial screw has been inserted in each main fragment, length and rotation are defined. Deformities can still be manipulated relatively well, whereas there are only limited correction possibilities for varus/valgus deformities. Therefore, it is recommended to insert the first screw in the distal fragment. The distal screws should be placed parallel to the knee joint. Then insert a screw in the distal fragment.

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Important: If a screw has to be removed and reinserted, use the torque-limiting screwdriver and not the power tool.

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Make a stab incision and insert the trocar through the drill sleeve. (Pic. 14)

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Determine screw length. The length of the condylar screws can be deduced from the table on page 05. Use screws of 26 mm length in the diaphyseal region.


Fig.14

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Options

In case of very thick cortex, pre-drill by using the pulling device or the drill bit 4.3mm.

The insertion of the initial screw tends to push the bone medially, especially in case of dense bone and/or unstable reductions.

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Determine screw length with Kirschner wire It is also possible to use the measuring device with a 2.0 mm Kirschner wire, placed through the centering sleeve (only use the gauge Ref.3420 with the guides Ref. 7016 or Ref.7017.Never use it with the guide Ref.7007 neither with Ref.7008). Using image intensification, insert the Kirschner wire to the desired depth leaving at least 5 mm between the tip of the Kirschner wire and the medial cortex. Measure the screw length over the Kirschner wire using the measuring device for Kirschner wires, leaving the centering sleeve in place, and round down to the nearest screw length. This will ensure that the tip of the screw will not protrude through the medial cortex (Pics. 15 and 16).


Fig.15

Fig.16

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Tighten the screw manually with the torque-limiting screwdriver Ref.1766. After one click, the optimum torque is reached. Insert a stopper into the Miplate insertion guide after screw insertion (Pics. 17 e 18)

Fig.17

Fig.18

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Note

To reduce the risk of stripping the screw head do not lock the screws at full speed. This can make it difficult to remove the implant. If the screw is difficult to insert or stops advancing prior to locking to the plate, remove the screw and clean the cutting flutes using a Kirschner wire. The screw can be reused if the socket has not been damaged. Should the screwdriver be difficult to remove after insertion, disconnect it from the power tool and remove the drill sleeve. After reconnecting the screwdriver to the power tool, withdraw the screwdriver from the screw.

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Predrill screw hole

Remove the trocar and thread the drill sleeve into the plate hole through the drill sleeve for Miplate insertion guide. Carefully drill the screw hole using the 4.3 mm drill bit (Pics. 19 e 20).

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Tighten the screw manually with the torque-limiting screwdriver Ref.1766. The optimum torque is reached after one click. Insert a stopper into the Miplate insertion guide after screw insertion to reduce the risk of stripping the screw head do not lock the screws at full speed. This can make it difficult to remove the implant. For long screws and thick cortical bone, ensure sufficient cooling during insertion.

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
Insert and lock the screw with the torque-limiting screwdriver through the drill sleeve for Miplate insertion guide (Pic. 18).

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Screwing the knurled nut onto the pulling device allows the bone to pull towards the plate. Since the tip of this instrument has a diameter of 4.0 mm, replacing it with a 5.0 mm locking screw still ensures good purchase in the bone (Pics.21and 22).

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Note:

It is important to monitor the advance of the screw tip carefully when inserting the pulling device. Stop the power tool before the pulling device is seated on the plate. Failure to do so may result in stripping the thread in the bone.

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If the reduction of the fracture causes difficulties, insert a Schanz screw antero-medially in the distal fragment, and use the screw as a joystick. The insertion of a Schanz screw into the proximal fragment can also be very beneficial. Should it still be impossible to perform a correct reduction, improve the access by enlarging the soft-tissue opening.

Fig.21


Fig.22

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Bending and twisting of the plate is not recommended as it may result in a misalignment between the holes of the insertion guide and the corresponding plate holes.

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Should the plate lie too ventral or too dorsal, the screws cannot be centred in the medullary canal. This position may compromise screw purchase

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Both screwdriver shaft and torque-limiting screwdriver are equipped with a self-holding mechanism. Apply slight pressure on pick-up to ensure that the screwdriver shaft penetrates the recess of the screw head.

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Should the screwdriver be difficult to remove after insertion, disconnect it from the power tool and remove the drill sleeve. After reconnecting the screwdriver to the power tool, withdraw the screwdriver from the screw.

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
Standard 4.5 mm cortex screws can be used through the insertion guide if required. Note that cortex screws cannot be inserted through the drill sleeve for Miplate insertion guide.

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If it is necessary to inset a screw in hole A, fix the insertion guide into the plate in an adjacent hole that was not previously used or if all the distal holes were occupied by a screw, pre drill hole A, remove the guide and thread the screw in a free hand way. Be careful with screw deviations to avoid damaging the cannal that was previously made by the drill.

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To ensure stability of the construct, the most proximal screw should be inserted last, just before removing the insertion guide. Remove the stabilization bolt and insert the screw through the drill sleeve. If hole A is unoccupied, it must be closed with a Screw Hole Insert to facilitate the application of the insertion guide for removing the implant.


Posicionamento Correto

Posicionamento Comprometido


Placa Distal de Fêmur Miplate

Miplate Femoral Distal Plate

PARAFUSO CORTICAL AUTO ROSQUEANTE Ø 4,5 mm
Self-Tapping Cortical Screw

REF. Titânio	REF. Inox		REF. Titânio	REF. Inox	
<i>Titanium</i>	<i>Stainless Steel</i>	Ø mm	<i>Titanium</i>	<i>Stainless Steel</i>	Ø mm
1502	6344	4,5 x 18	1511	6353	4,5 x 36
1503	6345	4,5 x 20	1512	6354	4,5 x 38
1504	6346	4,5 x 22	1513	6355	4,5 x 40
1505	6347	4,5 x 24	1514	6356	4,5 x 42
1506	6348	4,5 x 26	1515	6357	4,5 x 44
1507	6349	4,5 x 28	1516	6358	4,5 x 46
1508	6350	4,5 x 30	1517	6359	4,5 x 48
1509	6351	4,5 x 32	1518	6360	4,5 x 50
1510	6352	4,5 x 34			

PARAFUSO CORTICAL AUTO ROSQUEANTE Ø 5 mm
Self-Trapped Cortical Screw

REF. Titânio	REF. Inox		REF. Titânio	REF. Inox	
<i>Titanium</i>	<i>Stainless Steel</i>	Ø mm	<i>Titanium</i>	<i>Stainless Steel</i>	Ø mm
1473	6682	5,0 x 18	1486	6695	5,0 x 44
1474	6683	5,0 x 20	1487	6696	5,0 x 46
1475	6684	5,0 x 22	1488	6697	5,0 x 48
1476	6685	5,0 x 24	1489	6698	5,0 x 50
1477	6686	5,0 x 26	1490	6699	5,0 x 55
1478	6687	5,0 x 28	1491	6700	5,0 x 60
1479	6688	5,0 x 30	1492	6701	5,0 x 65
1480	6689	5,0 x 32	1493	6702	5,0 x 70
1481	6690	5,0 x 34	1494	6703	5,0 x 75
1482	6691	5,0 x 36	1495	6704	5,0 x 80
1483	6692	5,0 x 38	1496	6705	5,0 x 85
1484	6693	5,0 x 40	1497	6706	5,0 x 90
1485	6694	5,0 x 42			

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Placa Distal de Fêmur Miplate - Instrumentais

Miplate Femoral Distal Plate - Instruments

1766	TORQUÍMETRO 3,5 Nm <i>Torque Wrench</i>
2062	MEDIDOR DE PARAFUSO <i>Gauge for Screw</i>
2063	FRESA CORTICAL 4,5 mm <i>Cortical Reamer</i>
2072	BROCA AÇO INOX 3,2 mm <i>Stainless Steel Drill</i>
2074	GUIA DE BROCA P/ PLACA 3,2 mm <i>Drill Guide For Plate</i>
3420	MEDIDOR CANULADO 7,0 mm <i>Cannulated Gauge</i>
7004	CHAVE SEXTAVADA 3,5 mm-LONGA <i>Long Hexagonal Wrench</i>
7007	GUIA Ø 4,5 MM COM ROSCA <i>Guide with thread Ø 4,5 mm</i>
7008	GUIA Ø 2,5 MM COM ROSCA <i>Guide with thread Ø 2,5 mm</i>
7009	CHAVE CURVA <i>Curved Wrench</i>
7010	CHAVE EM T <i>T Wrench</i>
7011	FIO ROSQUEADO Ø 2,0 mm x 350 mm <i>Threaded Wire</i>

7012	PUNÇÃO Ø 7,4 MM <i>Puncture Ø 7,4 mm</i>
7013	GUIA CANULADO Ø 7,5 (CURTO) <i>Cannulated Guide Ø 7,5 (Short)</i>
7014	TRAVE POSICIONADORA <i>Positioning Bar</i>
7016	GUIA CANULADO Ø 2,5 com rosca <i>Cannulated Guide Ø 2,5 with thread</i>
7017	GUIA CANULADO Ø 4,5 com rosca <i>Cannulated Guide Ø 4,5 with thread</i>
7018	GUIA CANULADO Ø 7,5 (LONGO) <i>Cannulated Guide Ø 7,5 (Long)</i>
7021	CHAVE SEXTAVADA 3,5 MM <i>Hexagonal Wrench</i>
7022	CABO AXIAL COM ENGATE RÁPIDO <i>Axial Cable with Fast Assembly</i>
7023	BROCA Ø 4,3 MM X 90 MM <i>Drill</i>
7024	BROCA Ø 4,3 MM X 15 MM <i>Drill</i>
7026	TAMPÃO <i>Plug</i>
7048	CHAVE UNHA DE GATO <i>Wrench</i>
7049	GUIA DE REDUÇÃO <i>Reduction Guide</i>